

Ellman Assay

Markus Hoenicka

mhoen@farmr1.med.uth.tmc.edu

Version: \$Date: 1999/11/11 21:36:54 \$ \$Revision: 1.3 \$

Status: DRAFT

Objective

Determination of total sulfhydryl groups, protein-bound sulfhydryl groups, and free sulfhydryl groups in biological samples using DTNB (Ellman's reagent). Sensitivity 50 μ M to 1000 μ M.

Sedlak J, Lindsay RH (1968) Anal. Biochem. 25:192-205.

Chemicals

5,5'-Dithio-bis(2-nitrobenzoic acid) (DTNB); MW 396.3; Sigma D-8130

Ethylenediaminetetraacetate, tetrasodium salt, dihydrate (EDTA); MW 416.21; 99 %; Sigma ED4SS

Glutathione (GSH); MW 307.3; 99 %; reduced; Sigma G-4251

Methanol; p.A.; J.T. Baker 8045

Trichloroacetic acid (TCA); 72 %; unknown unknown

Triethanolamine-HCl (TEA); MW 185.7; Merck 1.08357

Trizma Base; MW 121.1; Sigma T-1503

Devices

Microplate reader

Multichannel pipettors

Repeaters

RIA-tube racks

Solutions

Stock solutions

Tris-HCl pH 8.2

24.22 g	Trizma Base	final: 30 mM
---------	-------------	--------------

8.32 g	EDTA	final: 3 mM
--------	------	-------------

Dissolve in 800 ml MilliQ. Adjust pH with 3 N HCl to 8.2. Add MilliQ to make 1000 ml. Store at room temperature.

Tris-HCl pH 8.9

48.44 g	Trizma Base	final: 262 mM
---------	-------------	---------------

8.32 g	EDTA	final: 13 mM
--------	------	--------------

Dissolve in 800 ml MilliQ. Adjust pH with 3 N HCl to 8.9. Add MilliQ to make 1000 ml. Store at room temperature.

TEA

1.86 g TEA 20 mM

Add MilliQ to make 500 ml. Store at 4°C.

TCA

13.9 ml TCA 10 %

Stir carefully into 70 ml MilliQ. Add MilliQ to make 100 ml. Store at room temperature.

DTNB29.7 mg DTNB final: 150 μ M

Dissolve in 25 ml Methanol. Store at 4°C.

Fresh solutions

Glutathione stock

31.1 mg GSH 20 mM

Add TEA to make 5 ml. Prepare fresh daily and keep on ice until used.

Procedures

GSH standardAdd 950 μ l TEA to 50 μ l of Glutathione stock to make a 1 mM GSH stock. Prepare the following concentration series:

Table 1. GSH standard concentration series

Conc. (μ M)	0	50	100	250	500	1000
GSH 1 mM (μ l)	0	25	50	125	250	500
TEA (μ l)	500	475	450	375	250	0

Total-SH

Add the following components to RIA tubes:

20 μ l Sample/Standard
 75 μ l Tris·HCl pH 8.2
 25 μ l DTNB
 400 μ l Methanol

Spin down samples at 3000 x g for 5 min at room temperature. Using a multichannel pipettor, transfer 3 x 90 μ l supernatant into a flat-bottom microplate. Measure extinction at 412 nm.**Free SH**

Use a V-bottom microplate to set up the following reaction:

50 μ l Sample/Standard
 2 x 25 μ l TCA

Spin down for 15 min at 1500 rpm and room temperature and keep the supernatant for later use.

Set up the following reaction in a flat-bottom microplate:

200 μ l	Tris-HCl pH 8.9
20 μ l	DTNB
50 μ l	Supernatant from the TCA precipitation

Measure at 412 nm.

Protein-bound SH

Calculate the difference of the Total-SH and the Free SH.